# DOPPELFARBSTAB- HAEMOMETER

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# Hemometer „Assistent“

## Directions for use and description.

The Hemometer „Assistent“ with two Glass Colour Standards is manufactured in accordance with the regulations of the German Society for Internal Medicine (G.I.M.) and certified as such.

### Directions for use:

- 1) Fill the graduated measuring tube up to the bottom-graduation line (mark 2) with n/10 Hydrochloric Acid.
- 2) Clean thoroughly finger tip or lobe of the ear with ether and take a drop of blood with a flamed needle.
- 3) Suck 20 cmm blood into the capillary pipette precisely up to the mark, wipe the pipette-point and blow the blood into the measuring tube. You will achieve a good mixture of the liquid by repeated suction and blowing. The mixture will be darkbrown and clear after about one minute. The moment of mixing blood and Hydrochloric Acid must be noted on the watch.
- 4) Add water by means of the water pipette and mix with the glass-stirrer until the colour of the solution matches the colour of the test rods.
- 5) Read the result by diffused day-light exactly three minutes after adding the blood to the Hydrochloric Acid.

**Precision:** The Hemometer „Assistent“ is equipped with two nonfading glass colour standards. A special manufacturing process guarantees very precise measuring results. Our Hemometers are subject to the control office of the German Society for Internal Medicine. Each directions for use — pamphlet bears the stamp of the above mentioned control-office.

### Hemometer Measuring Tubes (Spare) only with inscription

„Assistent“. In accordance with the regulations of the German Society for Internal Medicine (GIM) (German medical weekly magazine 76 of 1951, edition No. 35, page 1089) our Hemometer measuring tubes are marked „Assistent“. In order to achieve correct results it is of great importance to use measuring tubes marked „Assistent“ as the graduation of our „Assistent“ measuring tubes is calibrated in accordance with the colour of our glass-colourstandards. By using other measuring tubes wrong results will be obtained. The German Society for internal Medicine (GIM) points out in the mentioned weekly magazine that only by using original measuring tubes can correct results be assured.

### Measuring Tubes with graduation for gram-reading

Our hemometer measuring tubes are provided with a scala for reading the hemoglobin values expressed as grammes contained in 100 ml of blood since the German Society for Internal Medicine decided to discard the formerly also used arbitrary system of Hemometer units or percents. The main reason for dropping the old arbitrary systems of Hemometer Units or Percents and for adopting the exact method of measuring the absolute Hemoglobin-content in grams per 100 ml (g %) by the German Society for Internal Medicine (GIM) is that hemoglobin-values are interpreted differently in different countries, even in different parts of the same country due to climatic differences etc. The adoption of the formerly known „gram-scale“ method as the standard-method from now on by the German Society for Internal Medicine (GIM) will remove the insecurity which the older methods involved. In compliance with decision of this eminent body only such firms are allowed to mark their Hemometers „GIM“ which manufacture measuring tubes with the grammes-scale. This decision is proved in the German medical weekly magazin No. 35 page 1074/1951 (by Prof. Dr. med. Ludwig Heilmeyer, Freiburg, and Dr. med. Helmut Kilchling, Freiburg). For comparison between the absolute hemoglobin concentration and the former values see the following table containing besides the absolute hemoglobin values (in gramm percentage) the corresponding values on the basis 14,8g % = 100 HE (Hemometer Units), 16 g % = 100 HE and 17 g % = 100 HE.

The normal Hemoglobin value varies between 14 g % and 18 g % and depends largely on age, sex and surrounding conditions.

**TABLE**

to compare gramm reading with percentage reading:

Absolute HB- Concen- tration	Relative Hemoglobin values			Absolute HB- Concen- tration	Relative Hemoglobin values		
	100 HE 14,8g % (HE = Hemometer Units)	100 HE 16 g %	100 HE 17 g %		100 HE 14,8g %	100 HE 16 g %	100 HE 17 g %
4,0g %	27 HE	25 HE	24 HE				
4,5g %	30 HE	28 HE	26 HE				
etc.	etc.	etc.	etc.				

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---

## Hemoglobin Contents and Number of Erythrocytes

If the Hb contents is indicated in gram percent, then the old conception of a colour index originating from the relative Hb values, must be abandoned. Those values are replaced by the term  $Hb_E$  which clearly indicated the relation between Hb contents and the number of erythrocytes.

$$Hb_E = \frac{\text{Hb contents in 100 ccm blood}}{\text{Erythrocytes in 100 ccm blood}}$$

According to Bürker,  $Hb_E$  indicates the average Hb contents of the individual erythrocyte.

The determination of  $Hb_E$  is extremely simplified by the nomogramme given aside: The found values of hemoglobin and erythrocytes are connected with each other by means of a ruler or a strip of paper and the  $Hb_E$  value is read from the third vertical in prolongation of this connecting line.

The former colour-index 1,0 is equivalent to a  $Hb_E$  value of  $32 \times 10^{-12} \text{ g}$ .

The normal range of  $Hb_E$  is between 28 and  $32 \times 10^{-12} \text{ g}$ .

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